

Table of contents

\\University of Malta

<table border="1"> <tr> <td>UMRI_FHS -2025-00587 valentina</td> </tr> </table>	UMRI_FHS -2025-00587 valentina
UMRI_FHS -2025-00587 valentina	

<table border="1"> <tr> <td> <table border="1"> <tr> <td>New Exam</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>New Exam</td> </tr> </table>	New Exam
<table border="1"> <tr> <td>New Exam</td> </tr> </table>	New Exam	
New Exam		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>New Program 3</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>New Program 3</td> </tr> </table>	New Program 3
<table border="1"> <tr> <td>New Program 3</td> </tr> </table>	New Program 3	
New Program 3		

<table border="1"> <tr> <td> <table border="1"> <tr> <td> AAKnee_Scout_18ch pd_space_sag_iso T2Map_anatomical ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions </td> <td> spair SAG </td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td> AAKnee_Scout_18ch pd_space_sag_iso T2Map_anatomical ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions </td> <td> spair SAG </td> </tr> </table>	AAKnee_Scout_18ch pd_space_sag_iso T2Map_anatomical ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions	spair SAG
<table border="1"> <tr> <td> AAKnee_Scout_18ch pd_space_sag_iso T2Map_anatomical ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions </td> <td> spair SAG </td> </tr> </table>	AAKnee_Scout_18ch pd_space_sag_iso T2Map_anatomical ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions	spair SAG	
AAKnee_Scout_18ch pd_space_sag_iso T2Map_anatomical ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions	spair SAG		

\\University of Malta\UMRI_FHS -2025-00587 valentina\New Exam\New Program 3\AAKnee_Scout_18 ch
TA: 17 sec Coil Selection: Manual Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Resolution - Acceleration

Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Weak

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Prescan
Image Filter	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Slices per Slab	120
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FoV Read	450 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3.3 ms
TE	1.37 ms
Averages	1
Concatenations	1
AutoAlign	Knee
Coil Elements	18K

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Slices per Slab	120
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FoV Read	450 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3.3 ms
Multi-Slice Mode	Sequential
Series	Ascending
Concatenations	1

Contrast - Common

TR	3.3 ms
TE	1.37 ms
Flip Angle	12 deg
Fat-Water Contrast	Standard
Contrasts	1
Reconstruction	Magnitude

Geometry - AutoAlign

Slab Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	Knee
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Time to Center	7.5 s

Resolution - Common

FoV Read	450 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	224
Phase Resolution	100 %
Slice Resolution	68 %
Trajectory	Cartesian

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Manual
MSMA	S - C - T
Sagittal	Med >> Lat
Coronal	A >> P
Transversal	H >> F
Coil Combination	Adaptive Combine
Matrix Optimization	Off

Resolution - Acceleration

Acceleration mode	GRAPPA
Reference Scans	Integrated

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

Inline - MapIt

MapIt	None
Flip Angle	12 deg
Measurements	1
Contrasts	1
TE	1.37 ms
TR	3.3 ms
Save Original Images	On

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Bandwidth	540 Hz/Px
Asymmetric Echo	Weak

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

Sequence - Part 2

Introduction	On
RF Spoiling	On

System - Tx/Rx

Frequency 1H	123.243293 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Assistant

SAR Assistant	Off
---------------	-----

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Dynamic

Dynamic Mode	Standard
Flip Angle	12 deg
Measurements	1
Time to Center	7.5 s

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - Cardiac

Save Original Images	On
Contrasts	1
TE	1.37 ms
TR	3.3 ms

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Inline Composing	Off
------------------	-----

\\University of Malta\UMRI_FHS -2025-00587 valentina\New Exam\New Program 3\pd_space_sag_iso spair
TA: 6:38 min Coil Selection: Auto Voxel Size: 0.5×0.5×0.5 mm ³ Acc.: 4 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
Slices per Slab	240
Phase Oversampling	30 %
Slice Oversampling	33.3 %
FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	0.50 mm
TR	900.0 ms
TE	28.00 ms
Averages	1.0
Concatenations	1
AutoAlign	Knee > Standard
Coil Elements	18K

Contrast - Common

TR	900.0 ms
TE	28.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle Mode	Constant
Flip Angle	120 deg
Fat-Water Contrast	SPAIR
Fat Saturation	Strong
Dark Blood	Off
Blood Suppression	Off
Wrap-up Magn.	Restore
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement
Reordering	Radial

Resolution - Common

FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	0.50 mm
Base Resolution	320
Phase Resolution	100 %

Resolution - Common

Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration mode	CAIPIRINHA
Total Factor	4
Reference Scans	GRE/Separate
Acceleration Factor PE	2
Reference Lines PE	36
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Prescan
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
Slices per Slab	240
Phase Oversampling	30 %
Slice Oversampling	33.3 %
FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	0.50 mm
TR	900.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
AutoAlign	Knee > Standard
Initial Position	L1.4 A2.0 F8.0
L	1.4 mm
A	2.0 mm
F	8.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
-------------------	-----

Geometry - Tim Planning Suite

Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	ACS Restricted
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	H >> F
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Rotation	1.42 deg
A >> P	160 mm
F >> H	160 mm
R >> L	120 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.243293 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	900.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	SPAIR
Magn. Preparation	None
Dark Blood	Off
FoV Read	160 mm
FoV Phase	100.0 %
Phase Resolution	100 %
Dynamic Mode	Standard

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
----------	-----

Inline - Subtraction

Measurements	1
StdDev	Off
Save Original Images	On

Inline - Cardiac

Magn. Preparation	None
Save Original Images	On
TE	28.00 ms
TR	900.0 ms

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Inline Composing	Off
------------------	-----

Sequence - Part 1

Sequence Name	spcR
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Fast
Flow Compensation	None
Reordering	Radial
Bandwidth	422 Hz/Px
Echo Spacing	5.58 ms
Turbo Factor	60
Echo Train Duration	340 ms

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\University of Malta\UMRI_FHS -2025-00587 valentina\New Exam\New Program 3\T2Map_anatomica I SAG

TA: 4:45 min Coil Selection: Auto Voxel Size: 0.4×0.4×3.0 mm³ Acc.: 2 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	32
Distance Factor	10 %
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FoV Read	170 mm
FoV Phase	100.0 %
Slice Thickness	3.0 mm
TR	2440.0 ms
TE 1	13.80 ms
TE 2	27.60 ms
TE 3	41.40 ms
TE 4	55.20 ms
TE 5	69.00 ms
Averages	1
Concatenations	1
AutoAlign	Knee > Standard
Coil Elements	18K

Contrast - Common

TR	2440.0 ms
TE 1	13.80 ms
TE 2	27.60 ms
TE 3	41.40 ms
TE 4	55.20 ms
TE 5	69.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	180 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	5
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FoV Read	170 mm
FoV Phase	100.0 %

Resolution - Common

Slice Thickness	3.0 mm
Base Resolution	384
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	24
Phase Partial Fourier	4/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Off
Image Filter	Off

Geometry - Common

Slice Group	1
Slices	32
Distance Factor	10 %
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FoV Read	170 mm
FoV Phase	100.0 %
Slice Thickness	3.0 mm
TR	2440.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
AutoAlign	Knee > Standard
Initial Position	L1.4 A2.0 F8.0
L	1.4 mm
A	2.0 mm
F	8.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	ACS Restricted
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	H >> F
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	R69.0 A16.6 H43.7 mm
! Orientation	S > T7.4 > C6.0
! Rotation	1.42 deg
! A >> P	160 mm
! F >> H	160 mm
! R >> L	120 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
---------	----------

System - Tx/Rx

Frequency 1H	123.243293 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	2440.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FoV Read	170 mm
FoV Phase	100.0 %
Phase Resolution	100 %
Dynamic Mode	Standard

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - Cardiac

Magn. Preparation	None
Save Original Images	On

Inline - Cardiac

Contrasts	5
TE 1	13.80 ms
TE 2	27.60 ms
TE 3	41.40 ms
TE 4	55.20 ms
TE 5	69.00 ms
TR	2440.0 ms

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Inline Composing	Off
------------------	-----

Inline - MapIt

MapIt	T2 Map
Flip Angle	180 deg
Measurements	1
Contrasts	5
TE 1	13.80 ms
TE 2	27.60 ms
TE 3	41.40 ms
TE 4	55.20 ms
TE 5	69.00 ms
TR	2440.0 ms
Noise Threshold	15
Save Original Images	On

Sequence - Part 1

Sequence Name	se
RF Pulse Type	Low SAR
Gradient Mode	Fast
Bandwidth	228 Hz/Px

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\University of Malta\UMRI_FHS -2025-00587 valentina\New Exam\New Program 3\ep2d_diff_mddw_64_p2_s2 - PA ISOTROPIC sag

TA: 24 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc.: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	32
Distance Factor	30 %
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
TR	4500.0 ms
TE	54.00 ms
Concatenations	1
AutoAlign	Knee > Standard
Coil Elements	18K

Contrast - Common

TR	4500.0 ms
TE	54.00 ms
MTC	Off
Magn. Preparation	None
Fat-Water Contrast	Fat Saturation
Fat Saturation	Strong
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	80
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration mode	None
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
------------	-----

Resolution - Filter

Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	32
Distance Factor	30 %
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
TR	4500.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	P >> A
AutoAlign	Knee > Standard
Initial Position	L1.4 A2.0 F8.0
L	1.4 mm
A	2.0 mm
F	8.0 mm
Initial Orientation	Sagittal
Initial Rotation	-180.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	H >> F
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto

System - Adjustments

Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

Sequence - Part 2

Phase Correction	Internal
------------------	----------

System - Adjust Volume

Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Rotation	-178.58 deg
A >> P	160 mm
F >> H	160 mm
R >> L	83 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.243293 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	4500.0 ms
Concatenations	1

Physio - PACE

Resp. Control	Off
Concatenations	1

Diff

Diffusion Mode	MDDW
Diff. Directions	64
Diffusion Scheme	Monopolar
Diff. Weightings	1
b-value	0 s/mm ²
Averages	4
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Sequence - Part 1

Sequence Name	epse
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Performance
Bandwidth	2016 Hz/Px
Echo Spacing	0.60 ms
Free Echo Spacing	Off
Optimization	None
EPI Factor	80

Sequence - Part 2

Introduction	On
--------------	----

\\University of Malta\UMRI_FHS -2025-00587 valentina\New Exam\New Program 3\ep2d_diff_mddw_64_p2_s2 - AP ISOTROPIC sag 64 diff directions

TA: 5:12 min Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	32
Distance Factor	30 %
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
TR	4500.0 ms
TE	54.00 ms
Concatenations	1
AutoAlign	Knee > Standard
Coil Elements	18K

Contrast - Common

TR	4500.0 ms
TE	54.00 ms
MTC	Off
Magn. Preparation	None
Fat-Water Contrast	Fat Saturation
Fat Saturation	Strong
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	80
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration mode	None
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
------------	-----

Resolution - Filter

Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	32
Distance Factor	30 %
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FoV Read	160 mm
FoV Phase	100.0 %
Slice Thickness	2.0 mm
TR	4500.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Phase Encoding Dir.	A >> P
AutoAlign	Knee > Standard
Initial Position	L1.4 A2.0 F8.0
L	1.4 mm
A	2.0 mm
F	8.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	H >> F
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto

System - Adjustments

Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

Sequence - Part 2

Introduction	On
Phase Correction	Internal

System - Adjust Volume

Position	R69.0 A16.6 H43.7 mm
Orientation	S > T7.4 > C6.0
Rotation	1.42 deg
A >> P	160 mm
F >> H	160 mm
R >> L	83 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.243293 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	4500.0 ms
Concatenations	1

Physio - PACE

Resp. Control	Off
Concatenations	1

Diff

Diffusion Mode	MDDW
Diff. Directions	64
Diffusion Scheme	Monopolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	500 s/mm ²
Averages 1	4
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Sequence - Part 1

Sequence Name	epse
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Performance
Bandwidth	2016 Hz/Px
Echo Spacing	0.60 ms
Free Echo Spacing	Off
Optimization	None
EPI Factor	80